ANTITROMBINA III (AT III)

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ANTITHROMBIN III (AT III)

QUANTITATIVE DETERMINATION BY IMMUNOTURBIDIMETRY

For "in vitro" diagnostic

Summary

The primary function is to inactivate the proteolytic enzyme thrombin, prevent uncontrolled clot formation. AT III also functions to inhibit activated proteases such as coagulation factors IX, X, XI and XII.

Low levels of serum AT III are found in individuals suffering from severe hepatic disorders and as the result of protein-losing syndrome, particularly in the protein-losing nephropathies.

Diseases accompanied by inflammation may show elevated serum levels of Antithrombin III.

Principle

The reagent allows to quantify the level of human AT III present in the sample by comparing the turbidimetric response produced with that obtained from standard curve of known concentrations of human AT III. The reagent is an anti-serum anti-human AT III that reacts with the AT III of the serum giving protein aggregates. This aggregation process produces an increase in the Abs. of the system.

Reagents and Controls

Kit	Ref. 99 66 06
Contents	
1 x 2,8 ml anti-human AT III serum	Ref. 99 66 16
1 x 23,0 ml Buffer solution.	Ref. 99 66 46

Protein Calibrator.	Ref. 99 99 39
1 x 1 ml Pool of sera samples calibrated in Antithrombin III, IgA, IgG, IgM and Transferrin	
The concentration is stated on the enclosed sheet.	

Positive Control of Proteins	Ref. 99 66 86
1 x 1 ml sera sample positive for Antithrombin III, IgA, IgG, IgM and Transferrin	
The concentration is stated on the enclosed sheet.	

Caution

Human plasma used in the control preparation have been found negative in the reaction to HBsAg and HIV I/II. However they should always be handled with care.

All reagents contain sodium azide 0.09% as preservative. Waste products must be handled as per local regulations

Storage and stability

The antisera reagent, the buffer solution and the Protein Controls are stable up to the expiry date indicated on the label when kept at 2-8°C. Do not freeze antisera or buffer solution.

Sample

Recent plasma or kept at 2-8°C for not more than 48h. If the test should be performed later, it is recommended to freeze the serum. Avoid successive freezing and thawing. Discard hemolyzed or contaminated samples.

Quality Control

Calibrator or positive control should be included in each test series.

Procedure

To be used in clinical chemistry analyzers.

Use the antibody and the buffer solution directly. Shake gently the reagents before using it.

Use the sera samples without dilution.

Method

Volume of buffer solution: 280 µl
Volume of antibody solution: 35 µl
Volume of sample: 5 µl

Reaction time: 300 sec.
Reaction temperature: 37°C
Wavelength: 340 nm

Reaction range: up to 65 mg/dl

The before mentioned method is only indicative as the method should be specifically adapted to each autoanalyzer.

Calibration.

Prepare the standard curve using the protein calibrator. Use saline solution for concentration of 0 mg/dl.

Reference values

23 -39 mg/dl. It is recommended that each user calculates its own values.

Performance characteristics

The analytical characteristics depend on the instrument used. The mentioned results have been obtained using a TARGA analyzer. It is recommended to calculate these data for each instrument used.

VC inter serie: 2,1 %; VC intra series: 4,3 %

Prozone limit: No prozone effect was observed up to 300 mg/dl
Interferences: No interference was observed by Heparin (50 mg/dl), Sodium citrate (1000 mg/dl), Hemoglobin (1000 mg/dl), Bilirubin (20 mg/dl), Triglycerides (2500 mg/dl)

Limitations of the procedure

Samples that give an Abs. higher than that obtained for the last point in the standard curve should be diluted in saline (NaCl 0.9%) and the determination repeated.

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