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ANTICOAGULANTS. ANCILLARY PRODUCTS FOR COAGULATION TESTS.

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1 x 25mL Ref. 99 68 21

PRODUCT DESCRIPTION

Anticoagulant for use in hematology and biochemistry. Specially recommended for the determination of glycohemoglobin and for hematologic recounting. The blood collected with EDTA shows stability with cellular components and no signs of hemolysis for up to 8 days from collection. One drop is enough to inhibit coagulation of 5mL of total blood. The use of higher amounts of anticoagulant may yield falsely decreased hematocrit values.

REAGENT COMPOSITION

EDTA solution (0,738 mol/L) pH 8,4

STORAGE AND STABILITY

Store at Room Temperature. The reagent will remain stable until the expiration date stated on the label.

REAGENT PREPARATION

Reagent is ready to use.

CAUTION

Waste products must be handled as per local regulations.

SODIUM CITRATE (3,2%)

1 x 500mL Ref. 99 59 59

PRODUCT DESCRIPTION

Anticoagulant of choice for routine coagulation tests. Generally used for the determination of the VSG and the prothrombin time, the antithrombin III and the fibrinogen. The adequate blood/anticoagulant ratio is 9 parts of blood and 1 part of sodium citrate. Other citrate concentrations than the one specified will affect coagulation time.

REAGENT COMPOSITION

sodium Citrate solution 3,2%

STORAGE AND STABILITY

Store at 4-30°C. The reagent will remain stable until the expiration date stated on the label.

REAGENT PREPARATION

Reagent is ready to use.

CAUTION

Contains sodium azide (0.09%). Waste products must be handled as per local regulations.

CALCIUM CHLORIDE (0.025M)

1 x 100mL Ref. 99 88 02

PRODUCT DESCRIPTION

Product specially recommended for the Activated Partial Thromboplastin Time (APTT) and Prothrombin Time (PT) determinations.

REAGENT COMPOSITION

Calcium chloride solution 0,025 mol/L

STORAGE AND STABILITY

Store at Room Temperature. The reagent will remain stable until the expiration date stated on the label.

REAGENT PREPARATION

Reagent is ready to use.

CAUTION

Contains sodium azide (0.09%). Waste products must be handled as per local regulations.

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GB - EDTA K3

Danger

Hazard: Causes severe skin burns and eye damage.

Precautionary: Do not breathe dust/fume/gas/mist/vapours/spray. Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Specific treatment.

Contains: caustic potash,potassium hydroxide

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