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CLSI Guidelines and Standards, CLSI, Wayne, P.A

QUÍMICA CLÍNICA APLICADA S.A.

Empresa Certificada ISO 9001 / ISO 13485

A7 Km 1081 - P.O. Box 20 - E43870 AMPOSTA / SPAIN

Tel. ++ 34 (977) 70. 62. 30 Fax ++ 34 (977) 70. 30. 40

Revisión: 05.2018

PRO4-9_COLCITO_5

STAINS FOR HISTOCYTOLOGY. CITOPLASM STAINING

For in vitro diagnostic

PRINCIPLE

The Papanicolaou test, which owes its name to Dr Georgios Papanicolaou, a pioneer in cytology and cancer, is performed to diagnose cervical cancer or cancer in other organs. It makes it possible to detect changes in cells which may be cancer precursors. The first step in this technique is the staining of nuclei with haematoxylin. The second step is the staining of cytoplasm with an orange solution, which stains the mature cells and the keratinised cells with different intensities. The so-called polychrome solution, which is a mixture of eosin, light green SF and Bismarck brown, is used in the third step. The differentiation of simple squamous epithelium is shown with the polychrome solution. The Haematoxylin and Eosin stain is the most commonly used stain in histology.

The dyes for the staining of the cytoplasm are fixed in the extracellular components of the tissue and enable its differentiation. Eosin is the main component of this group. The stain spreads easily between the tissue structures and is fixed, due to its acidic nature, to the basic radicals present in the tissue proteins. The Orange G stain is also an acid stain which has greater penetration capacity than eosin in dense structures, therefore enabling a greater staining of keratinised cells. Stains from the EA group combine the staining capacity of an acidic dye, such as Eosin, with the staining obtained using the Light Green stain, which, due to its less acidic nature, adds staining shades to the cytoplasmic structure.

DIAGNOSTIC USE

For cytoplasmic staining in histocytology.

REAGENTS

Alcoholic eosin 1%	1 x 1000 mL	Ref. 99 14 15
Composition		
Eosin	1%	
Eosin 1%	1 x 1000 mL	Ref. 99 30 06
Composition		
Eosin aqueous solution	1%	
Papanicolaou OG-6	1 x 500 mL 1 x 1000 mL 1 x 2500 mL 1 x 5000 mL	Ref. 99 46 45 Ref. 99 70 10 Ref. 99 76 02 Ref. 99 74 20
Composition		
Orange G	0.25%	
Ethanol	90%	
Papanicolaou EA-36	1 x 1000 mL 1 x 5000 mL	Ref. 99 69 89 Ref. 99 55 36
Composition		
Eosin	0.20%	
Light Green	0.1%	
Methanol	90%	
Papanicolaou EA-50:	1 x 500 mL 1 x 1000 mL 1 x 2500 mL 1 x 5000 mL	Ref. 99 19 07 Ref. 99 70 12 Ref. 99 83 34 Ref. 99 19 22
Composition		
Eosin	0.20%	
Light Green	0.05%	
Ethanol	80%	
Methanol	10%	
Papanicolaou EA-65:	1 x 1000 mL	Ref. 99 69 64
Composition		
Eosin	0.15%	
Light Green	0.02%	
Methanol	90%	

PREPARATION OF THE WORKING REAGENT:

The stain is ready for use.

STORAGE AND STABILITY

Reagents which are stored at 15-30°C and protected from light are stable until the expiry date stated on the label. Containers must always be kept tightly closed.

Over time, a light precipitate may form in some reagents. This does not affect their functionality. The stain should be filtered before use.

MATERIAL REQUIRED (NOT SUPPLIED)

General-purpose pathological anatomy laboratory material.

Haematoxylin stain.

Different concentrations of ethanol. Xylene or Xylene Substitute. Mounting medium

Microscope

PRECAUTIONS

The products are for professional use. All the reagents must be handled with care by trained technicians. The safety data sheet should be consulted before using the reagents.

Waste disposal must be carried out according to the local regulations in force.

SAMPLE

Cytological samples of different origins: gynaecological, sputum, liquid biopsies.

Histological samples.

Handling of the samples must be carried out in accordance with the established protocols in each laboratory for the preparation of samples for staining with the Papanicolaou method or with the Haematoxylin and Eosin stain.

Handle the samples with care due to their potentially infectious nature.

QUALITY CONTROL

Following the quality control practices defined by the CLSI (formerly NCCLS) is recommended.

Suitable controls must be carried out for each set of stains to avoid inaccurate results.

PROCEDURE

PAPANICOLAOU STAIN

1. The smears which have been spray-fixed should be immersed in water for 5-10 minutes so that the carboxymethyl cellulose (CMC) that they contain is completely dissolved
2. Immerse in Harris's Haematoxylin for 1 minute
3. Rinse with water for 15 seconds
4. Rinse with ethanol-HCl (0.25%) for 10 seconds
5. Rinse with water for 15 seconds
6. Rinse with ammonia water (0.05% in NH₃) for 5 seconds
7. Rinse with water for 15 seconds
8. Immerse in ethanol 96° for 10 seconds (twice, using a different container each time)
9. Immerse in OG-6 stain for 3 minutes
10. Rinse in ethanol 96° for 15 seconds (twice, using a different container each time)
11. Immerse in EA 50, EA 36 or EA 65 stain for 3 minutes
12. Rinse in ethanol 96° for 15 seconds (twice, using a different container each time)
13. Immerse in absolute ethanol for 15 seconds (twice, using a different container each time)
14. Immerse in ethanol/xylol (50%) for 15 seconds
15. Immerse in xylol for 15 seconds
16. Immerse in xylol for 2 minutes. Mount the preparation

HAEMATOXYLIN AND EOSIN STAIN

1. Dewax the preparation
2. Rinse with distilled water for 1 minute after the last immersion in ethanol
3. Immerse in Haematoxylin for 8 minutes
4. Rinse with water for 1 minute
5. Rinse with ethanol-HCl (0.5%) for 13 seconds
6. Rinse with water for 15 seconds
7. Rinse with ammonia water (0.05% in NH₃) for 30 seconds
8. Rinse with water for 1 minute
9. Immerse in Eosin 1% for 30 seconds
10. Immerse in ethanol 96° for 6 seconds (twice, using a different container each time)
11. Immerse in absolute ethanol for 6 seconds (three times, using a different container each time)
12. Immerse in xylol/eucalyptol for 6 seconds (twice, using a different container each time)
13. Mount the preparation

RESULTS

Nuclei: blue-violet

Acidophilic cytoplasm: reddish

Basophilic cytoplasm: blue

Keratinised cytoplasm: orange

NOTES

The staining intensity is proportional to the staining time. The stain must be kept away from moisture in order to obtain optimal results.

Each user may apply the different versions of this procedure, both manual and automated, adapting it to their standard method.

REFERENCES

Marshall, P.N., Galbraith, W., Bacus, J.W., (1979). Anal. Quant. Cytology, 1, 169-178.

Baker, J.R., (1962). Quart. J. Micr. Sci., 103, 493-517.

Boon, M.E., Drijver, J.S., Routine Cytological staining techniques, 1st. ed (1986)

CLSI Guidelines and Standards, CLSI, Wayne, P.A

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Révision: 05.2018

PRO4_9_COLCITO_5

EOSINA ALCOHOLICA 1% ALCOHOLIC EOSIN 1%

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GB - EOSIN ALCOHOLIC 1%

Danger

Hazard: Highly flammable liquid and vapour.

Precautionary: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. In case of fire: Use the means described in point 5 of the Safety Data Sheet. Store in a well-ventilated place. Keep cool. Dispose the contents/container according to point 13 of the Safety Data Sheet.

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GB - PAPANICOLAOU (EA-36)

Danger

Hazard: Highly flammable liquid and vapour. Toxic if swallowed, in contact with skin or if inhaled. Causes damage to organs.

Precautionary: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not breathe dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF exposed or concerned: Call a POISON CENTER/doctor/... Specific treatment. In case of fire: Use the means described in point 5 of the Safety Data Sheet. Store in a well-ventilated place. Keep container tightly closed.

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GB - PAPANICOLAOU (EA-50)

Danger

Hazard: Highly flammable liquid and vapour. Harmful if swallowed, in contact with skin or if inhaled. Causes damage to organs.

Precautionary: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Do not breathe dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF exposed or concerned: Call a POISON CENTER/doctor/... Specific treatment. In case of fire: Use the means described in point 5 of the Safety Data Sheet.

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Contient: méthanol

PAPANICOLAOU EA-65 POLYCHROMIC SOLUTION (EA-65)

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