### PRINCIPIO

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Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

# FIBRINOGEN

## CLAUSS METHOD

For "in vitro" diagnostic

### PRINCIPLE

The fibrinogen determination by clotting time is based on the method originally described by Clauss. In the presence of an excess of thrombin, the fibrinogen is transformed into fibrin and the time to clot formation is inversely proportional to the concentration of fibrinogen present in the plasma sample.

The QCA Thrombin reagent placed in contact with the plasma of the patient, acts on the fibrinogen to be converted into fibrin. This fibrin polymerizes to form a dense network which contributes to the clot formation. The training time can be measured manually or automatically.

### DIAGNOSTIC USE

Fibrinogen is a plasma glycoprotein which is involved in the coagulation processes. Fibrinogen increases in inflammatory processes and tissue diseases (infarction, leukaemia and lupus). Fibrinogen decreases in severe hepatic diseases. Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

### REAGENTS

#### KIT Ref. 99 94 30

Content:	
A. 4 x 2 mL Thrombin	Ref. 99 94 10
B. 1 x 50 mL Buffer	Ref. 99 94 11
C. 1 x 1 mL Fibrinogen calibrator	Ref. 99 94 12

### WORKING REAGENT PREPARATION

Reagent A: rehydrate one vial with 2 mL of distilled water. Mix gently by inversion and wait 10 min for complete homogenization. Once open and rehydrated is stable 7 days at 2-8°C.  
Reagent B: imidazol buffer is ready to use.  
Reagent C: rehydrate one vial with 1 mL of distilled water. Mix gently by inversion and wait 20 min for complete homogenization. Once open and rehydrated is stable 8h at room temperature ( $\leq 25^{\circ}\text{C}$ ) or 2 days at  $-20^{\circ}\text{C}$ . The concentration is stated on the vial label.

### WORKING REAGENT COMPOSITION

A: Freeze-dried bovine thrombin buffered at pH 7.4 with preservatives and stabilizers.  
B: Imidazole buffer solution 15mM at pH 7.4, which contains NaCl 125mM with preservatives and stabilizers  
C: Freeze-dried pool of human plasma, standardized and assayed for fibrinogen.

### STORAGE AND STABILITY

When stored at 2-8°C, reagents are stable until the expiration date stated on the label. Do not use reagents over the expiration date. Do not freeze.

### Signs of reagent deterioration:

Presence of particles in the reagent B. Quality control values outside acceptance range.

### ADDITIONAL EQUIPMENT

General laboratory equipment.  
Coagulometer or bath 37°C and stopwatch.

### SAMPLE

Mix carefully 9 parts of fresh blood with 1 part of trihydrated trisodium citrate 3.2% (0.109M, Ref. 99 59 59). Centrifuge at 3.000 rpm/5 min and separate the supernatant plasma. Store at 2-8°C until prior to assay. (Not more than 2 hours)

### CAUTION

The safety statements are on the label. It is recommended to read SDS before reagent manipulation.  
The control plasma and patient plasma must be handled with care.  
Waste products must be handled as per local regulations.

### QUALITY CONTROL

Control plasma should be included in each test series.

#### Set of controls Ref. 99 25 50. Contents:

1 x 1 mL Coagulation Control, Low Level  
1 x 1 mL Coagulation Control, High Level  
Pool of human plasma with standardized levels of coagulation factors.

Each particular laboratory should establish its own control program and correctives actions of measure deviations.

### PROCEDURE

#### 1. Plasma dilution

Dilute plasma and control plasma 1:10 with buffer (reagent B) (1 plasma + 9 buffer)

#### 2. Method

- Preincubate thrombin (reagent A) at 37°C/10-15min
- Dispense 0.2 mL of diluted plasma into a test tube and incubate at 37°C during 2min
- Add 0.1 mL of thrombin (reagent A) at 37°C and start the stopwatch.
- Record the time for the clot to be formed.

#### 3. Calibration curve

- Rehydrate fibrinogen calibrator vial as described above.
- Prepare dilution series as following:

	150%	100%	50%	25%
Buffer (mL)	0.85	0.90	1.90	3.90
Calibrator (mL)	0.15	0.10	0.10	0.10

- Determine the coagulation time of each sample dilution according to the described method.

- Coagulation time obtained can be plotted on a graph vs each concentration in the calibration curve. Calibration curve allows to change into fibrinogen concentration (%) the plasma coagulation times of the patients samples by interpolation.

Calibration curve must be repeated at every fibrinogen lot change.

### CALCULATION

Each particular laboratory should establish its own calibration curve as described above, and then determine by interpolation the sample concentration. Duplicate coagulation time determinations are recommended, and then calculate the mean value of double determinations of each plasma sample before to interpolate in the calibration curve.

### REFERENCES VALUES

Each particular laboratory should establish its own normal range, obtained from samples of a representative population, using its own instrumentation, blood collection methods and assaying procedures.

It is advisable not to establish the mean value with less than 20 individuals. However and as a guideline, the values of a certain ambulatory population are related herewith:  
200-400 mg/dL (2-4 g/L)

### PERFORMANCE CHARACTERISTICS

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using mechanical coagulometer.

Sensitivity, as minimum detectable concentration: 5mg/dL

Accuracy: 99.5%

Repetitivity, as CV%: 5.25%

Reproducibility, as CV%: 5.9%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

### INTERFERENCES

Important interferences have been observed in samples with degradation fibrinogen products (FDP or D-Dimer).

Acute inflammatory reactions can increase circulating fibrinogen (factor I).

Hemolysis can cause clotting factor activation and end point detection interference.

Drugs and substances consumption may interfere (See references)

### REFERENCES

Clauss, A. Acta Haemat., 1957, 17, 237-246.

NCCLS, H30-A2, 2nd Ed., 2001 21(18).

Kitchen, S.; McCraw, A.; Echenagucia, M.; Diagnosis of Hemophilia and Other Bleeding Disorders, 2nd Ed.; 2010, 55-57.

Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

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