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For in vitro diagnostic

## PRINCIPLE

The Papanicolaou test, which owes its name to Dr Georgios Papanicolaou, a pioneer in cytology and cancer, is performed to diagnose cervical cancer or cancer in other organs. It makes it possible to detect changes in cells which may be cancer precursors. The first step in this technique is the staining of nuclei with haematoxylin. The second step is the staining of cytoplasm with an orange solution, which stains the mature cells and the keratinised cells with different intensities. The so-called polychrome solution, which is a mixture of eosin, light green SF and Bismarck brown, is used in the third step. The differentiation of simple squamous epithelium is shown with the polychrome solution. The Haematoxylin and Eosin stain is the most commonly used stain in histology.

Haematoxylin stains are basically formed from metallic complexes of oxidised haematoxylin with divalent and trivalent ions. The mechanism of staining is carried out between the covalent bonds of the metal-haematoxylin compound and the anionic radicals of the tissue. The anionic groups of phosphoric acid in DNA will intervene if there is nuclear material in the reaction. The selectivity of the stain increases with the acid medium of the dye, since it hinders the staining of amphoteric elements, such as the amino acids of proteins. The nuclei appear blue to dark violet.

Progressive and regressive staining can be performed with haematoxylin. The structures of the nuclei show greater differentiation in regressive staining and are more easily seen.

## DIAGNOSTIC USE

For the staining of nuclei in histocytology.

## REAGENTS

### Harris's Haematoxylin

	1 x 500 mL	Ref. 99 22 32
	1 x 1,000 mL	Ref. 99 59 60
	1 x 5,000 mL	Ref. 99 12 97

### Composition

Haematoxylin C.I. No. 75 290	0.5%
Aluminium Sulphate	6%
Ethylene glycol	15%

### Carazzi's Haematoxylin

	1 x 1,000 mL	Ref. 99 16 43
	1 x 5,000 mL	Ref. 99 41 75

### Composition

Haematoxylin C.I. No. 75 290	0.1%
Aluminium/Potassium Sulphate	6%
Glycerin	20%

### Mayer's Haematoxylin

	1 x 1,000 mL	Ref. 99 77 50
	1 x 5,000 mL	Ref. 99 06 38

### Composition

Haematoxylin C.I. No. 75 290	0.5%
Aluminium Ammonium Sulphate	5%
Glycerin	30%

## PREPARATION OF THE WORKING REAGENT:

The stain is ready for use.

## STORAGE AND STABILITY

Reagents which are stored at 15-30°C and protected from light are stable until the expiry date stated on the label. Containers must always be kept tightly closed.

Over time, a light precipitate may form in some reagents. This does not affect their functionality. The stain should be filtered before use.

## MATERIAL REQUIRED (NOT SUPPLIED)

General-purpose pathological anatomy laboratory material.

Dye for staining cytoplasm: EA 50, EA 36 or EA 65 stains, Orange G-6 stain, Alcoholic or aqueous solution of eosin at 1%

Different concentrations of ethanol.

Xylene or Xylene Substitute

Mounting medium

Microscope

## PRECAUTIONS

The products are for professional use. All the reagents must be handled with care by trained technicians. The safety data sheet should be consulted before using the reagents.

Waste disposal must be carried out according to the local regulations in force.

## SAMPLE

Cytological samples of different origins: gynaecological, sputum, liquid biopsies.

Histological samples.

Handling of the samples must be carried out in accordance with the established protocols in each laboratory for the preparation of samples for staining with the Papanicolaou method or with the Haematoxylin and Eosin stain.

Handle the samples with care due to their potentially infectious nature.

## QUALITY CONTROL

Following the quality control practices defined by the CLSI (formerly NCCLS) is recommended.

Suitable controls must be carried out for each set of stains to avoid inaccurate results.

## PROCEDURE

### PAPANICOLAOU STAIN

1. The smears which have been spray-fixed should be immersed in water for 5-10 minutes so that the carboxymethyl cellulose (CMC) that they contain is completely dissolved
2. Immerse in Harris's Haematoxylin for 1 minute
3. Rinse with water for 15 seconds
4. Rinse with ethanol-HCl (0.25%) for 10 seconds
5. Rinse with water for 15 seconds
6. Rinse with ammonia water (0.05% in NH<sub>3</sub>) for 5 seconds
7. Rinse with water for 15 seconds
8. Immerse in ethanol 96° for 10 seconds (twice, using a different container each time)
9. Immerse in OG-6 stain for 3 minutes
10. Rinse in ethanol 96° for 15 seconds (twice, using a different container each time)
11. Immerse in EA 50, EA 36 or EA 65 stain for 3 minutes
12. Rinse in ethanol 96° for 15 seconds (twice, using a different container each time)
13. Immerse in absolute ethanol for 15 seconds (twice, using a different container each time)
14. Immerse in ethanol/xylol (50%) for 15 seconds
15. Immerse in xylol for 15 seconds
16. Immerse in xylol for 2 minutes. Mount the preparation

### HAEMATOXYLIN AND EOSIN STAIN

1. Dewax the preparation
2. Rinse with distilled water for 1 minute after the last immersion in ethanol
3. Immerse in Haematoxylin for 8 minutes
4. Rinse with water for 1 minute
5. Rinse with ethanol-HCl (0.5%) for 13 seconds
6. Rinse with water for 15 seconds
7. Rinse with ammonia water (0.05% in NH<sub>3</sub>) for 30 seconds
8. Rinse with water for 1 minute
9. Immerse in Eosin 1% for 30 seconds
10. Immerse in ethanol 96° for 6 seconds (twice, using a different container each time)
11. Immerse in absolute ethanol for 6 seconds (three times, using a different container each time)
12. Immerse in xylol/eucalyptol for 6 seconds (twice, using a different container each time)
13. Mount the preparation

### MAYER'S HAEMATOXYLIN AND EOSIN STAIN

1. Dewax the preparation
2. Rinse with distilled water for 5 minutes after the last immersion in ethanol
3. Immerse in Mayer's Haematoxylin for 10 minutes
4. Rinse with water for 5 minutes
5. Immerse in alcoholic Eosin 1% for 2 minutes
6. Immerse in ethanol 96° for 10-15 seconds (twice, using a different container each time)
7. Immerse in absolute ethanol for 10-15 seconds (twice, using a different container each time)
8. Immerse in xylol (three times, using a different container each time)
9. Mount the preparation

## RESULTS

Nuclei: blue-violet  
Acidophilic cytoplasm: reddish  
Basophilic cytoplasm: blue  
Keratinised cytoplasm: orange

## NOTES

The staining intensity is proportional to the staining time. The stain must be kept away from moisture in order to obtain optimal results. Each user may apply the different versions of this procedure, both manual and automated, adapting it to their standard method.

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