#### PRINCIPIO

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Ratnoff, O. D., Crum, S. D. (1964). J. Lab. Clin. Med., 63, 359.

Brandt, J.T., Triplett, D.A. (1981). Amer. J. Clin. Path., 76, 530.

Procedure for the handling and processing of blood specimens, H18-A3, NCCLS, 24(38), 2004.

One-Stage Prothrombin Time Test and Activated Partial Thromboplastin Time Test, H47-A2, CLSI, 28(20), 2008.

**PRINCIPLE**

The reagent is formed by a phospholipid obtained from a chloroform extract of rabbit brain, that performs as a platelet substitute and ellagic acid which is used as the activator of the coagulation mechanism. Once optimal activation is achieved by incubation of plasma with the reagent (at 37°C/3min), the mixture was recalcified and the coagulation time is determined.

**DIAGNOSTIC USE**

The determination of Activated Partial Thromboplastin Time (APTT) is mainly used to detect deficiencies in the first stage of the coagulation mechanism, principally factors VIII, IX, XI, XII and Prekallikrein (Fletcher Factor). It is also sensitive to deficiencies of the rest of factors except Factor VII.

It is commonly used for monitoring heparin anticoagulant therapy.

Levels approximately 40% of the normal value of Factors VIII, IX, XI and XII are sufficient to produce normal values of APTT. For levels lower than 40%, the APTT will be longer. Heparin, in the presence of adequate level of Antithrombin III, will also prolong the clotting time value.

Single test result could not be used to make a clinical diagnosis. It should integrate clinical and laboratory data.

**REAGENTS**

**1 x 4 mL (Ref. 99 05 89)** Contents:  
 1 x 4 mL Reagent for APTT determination Ref. 99 01 87

**1 x 10 mL (Ref. 99 36 34)** Contents:  
 1 x 10 mL Reagent for APTT determination Ref. 99 02 87

**Kit 4 x 10 mL (Ref. 99 26 89)** Contents:  
 4 x 10 mL Reagent for APTT determination Ref. 99 02 87

**Kit 10 x 4 mL (Ref. 99 40 30)** Contents:  
 10 x 4 mL Reagent for APTT determination Ref. 99 01 87

**WORKING REAGENT PREPARATION**

Reagent is ready to use.

**REAGENT COMPOSITION**

Extract of rabbit phospholipid with ellagic acid as activator, stabilizers and preservatives.

**STORAGE AND STABILITY**

When stored at 2-8°C, the reagent is stable until the expiration date stated on the label if the bottle is properly closed and preserved from contamination.

**After a long period of storage, yellow sediment may appear. In such a case, mix by gentle inversion to assure a proper activator suspension.**

Do not use reagents over the expiration date. Do not freeze.

**Signs of reagent deterioration:**

Presence of particles in the reagent after suspension. Quality control values outside acceptance range.

**ADDITIONAL EQUIPMENT**

General laboratory equipment.  
 Coagulometer or bath 37°C and stopwatch.

**Calcium Chloride 0.025M 1 x 100 mL Ref. 99 88 02**

**SAMPLE**

Mix carefully 9 parts of fresh blood with 1 part of trihydrated trisodium citrate 3.2% (0.109M, Ref. 99 59 59). Centrifuge at 3.000 rpm/5min. and separate the supernatant plasma. Store at 2-8°C until prior to assay. (Not more than 2 hours).

**QUALITY CONTROL**

Control plasma should be included in each test series.

**Set of controls (Ref. 99 25 50).** Contents:

1 x 1 mL Coagulation Control, Low Level

1 x 1 mL Coagulation Control, High Level

Pool of human plasma freeze-dried with standardized levels of coagulation factors.

Each particular laboratory should establish its own control program and correctives action of measure deviations.

**PROCEDURE**

1. Preincubate the CaCl<sub>2</sub> 0.025M solution at 37°C.
2. Dispense 0.1 mL of plasma into a test tube.
3. Mix, by gentle inversion, the Hemoscann® reagent and add 0.1 mL to the plasma.
4. Incubate the mixture plasma-reagent during 3 min. (time of activation) at 37°C.
5. After this period of time, add 0.1 mL of the CaCl<sub>2</sub> solution at 37°C and start the stopwatch.
6. Record the time for the clot to be formed.

**RESULTS**

The results shall be expressed as clotting time in seconds. Calculate the mean value of duplicate samples and controls.

**REFERENCES VALUES**

Each particular laboratory should establish its own normal range, obtained from samples of a representative population, using its own instrumentation, blood collection methods and assaying procedures.

It is recommended a minimum of 20 individual samples to establish the mean value. The values obtained in samples from healthy individuals are indicated as a guideline

Manual	20 - 35s
Mechanical	22 - 36s

**PERFORMANCE CHARACTERISTICS**

The performance characteristics depend on the method used. It is recommended to calculate these data for each particular test protocol. These results have been obtained using mechanical coagulometer.

Sensitivity to Heparin:

Heparin concentration U/mL	0.0	0.1	0.2	0.3	0.4	0.6	0.7	0.8
Coagulation time sec.	24.6	27.7	33.5	39.6	49.4	68.9	88.5	108

Repetitivity, as CV%: 1.8%

Reproducibility, as CV%: 3.4%

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagent.

Details of the performance studies are available on request.

**INTERFERENCES**

Turbid, icteric, lipemic or grossly haemolysed samples will give erroneous results. APTT values can be shortened as a result of the response to an acute inflammatory reaction in which the elevation of Factor I (Fibrinogen) is very common.

Contraceptives, estrogens, pregnancy, coumarin derivatives, heparin, asparaginase and naloxone, have been reported to influence the APTT values.

**CAUTION**

The safety statements are on the label. It is recommended to read MSDS before reagent manipulation.

The control plasma and patient plasma must be handled with care.

Waste products must be handled as per local regulations.

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